



Mammography
Reporting System

IHE Integration Statement			
Vendor	Product Name	Version	Date
MRS, Inc.	MRS	6.5	31 January 2006
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Reporting Workflow (RWF)	Report Creator	none	
Patient Synchronized Applications	Patient Context Participant	none	
<u>Internet address for MRS IHE information:</u> www.mrsys.com/ihe			
Links to Standards Conformance Statements of the Implementation			
HL7			
DICOM	www.mrsys.com/Integration/DICOM_MW_SCU_Conformance_Statement.pdf www.mrsys.com/Integration/DICOM_SR_SCU_Conformance_Statement.pdf		
CCOW	www.mrsys.com/Integration/CCOW_Conformance_Statement.pdf www.mrsys.com/Integration/HL7_Context_Management_Standard.pdf		
Links to general information on IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	